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# integrity™ V500

Generación 2 (G2)

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- Ambient temperature of -20 °C to +50 °C (for 30 days)
- Ambient temperature of -20 °C to +40 °C (for 90 days)
- Ambient temperature of -20 °C to +30 °C (for 1 year)
- Relative humidity of 10% to 100%
- Atmospheric pressure range of 500 hPa to 1060 hPa

#### Integrity (sin Paquete de baterías)

- Ambient temperature of -40 °C to +70 °C
- Relative humidity of 10% to 100%
- Atmospheric pressure range of 500 hPa to 1060 hPa

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- Ambient temperature of 10 °C to 39 °C
- Relative humidity of 30 % to 75 %
- Atmospheric pressure range of 700 hPa to 1060 hPa

## PRECAUCIÓN

Following the "Guidance for Industry, FDA Reviewers and Compliance on Off-The-Shelf Software Use in Medical Device", issued by FDA on September 9, 1999, the customers should not use any software other than that specified, as it will violate the safety, effectiveness and design controls of our Integrity V500 and that such use may result in an increased risk to users and patients.

Unauthorized alteration of software provided, including downloading non-validated or unauthorized off-the-shelf software may void or breach existing service agreements and warranties. For additional information, contact Vivosonic Inc.'s Customer Support.

## PRECAUCIÓN

Please follow these recommendations to ensure the test data is accurate.

- Substitution of any supplied components with different components may result in measurement error.
- No other software may be installed onto the computer on which the Integrity™ software is installed.
- Handle the system with care.
- The system requires calibration of the transducers – A81/A82 Amplitrodes, ER-3A-800 Earphones, B71W Bone Conductor, OAE probes, and H-800 supra-aural headphones – annually on each anniversary from the date of manufacture, and every time a transducer is dropped, subjected to mechanical shock, or immersed in a liquid substance. Otherwise, stimuli presented to the patient as specified may lead to incorrect test results and misdiagnoses.
- DO NOT force the ear tip into the ear canal.
- Always carry this device with you when traveling to avoid mishandling and damage of stored or checked luggage.
- DO NOT ship this device for service in packaging other than the packaging supplied with the system from the manufacturer, or comparable packaging.
- Signal inputs and outputs are intended for analysis only in connection to the specified equipment described in this manual.
- The product DOES NOT require sterilization. Alcohol can be used to clean the transducers.
- This equipment does not have protection against water penetration (IPX0).
- The system can only be repaired by Vivosonic Inc. or authorized dealers of Vivosonic Inc.

## Immunity

The Integrity System is intended for use in the electromagnetic environment specified below. The customer or user of the Integrity should ensure that it is used in such an environment.

### Equipment and systems that are NOT life-supporting

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment- Guidance
			Portable and mobile communications equipment should be separated from the Integrity by no less than the distances calculated, listed below.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	V1=3 Vrms	$D=(3.5/V1)(\text{Sqrt } P)$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 MHz	E1=3 V/m	$D=(3.5/V1)(\text{Sqrt } P)$ 80 to 800 MHz
			$D=(7/E1)(\text{Sqrt } P)$ 800 MHz to 2.5 GHz Where P is the max. power (watts) and D is the recommended separation distance (meters). Field strength from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1). Interference may occur in the vicinity of equipment containing a transmitter.

### Immunity for all equipment and systems

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment- Guidance
ESD IEC 61000-4-2	±6kV Contact ±8kV Air	±6kV Contact ±8kV Air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the r/h should be at least 30%.
EFT IEC 61000-4-4	±2kV Mains ±1kV I/Os	N/A. Battery Operated Device	N/A
Surge IEC 61000-4-5	±1kV Differential ±2kV Common	N/A. Battery Operated Device	N/A
Voltage Dips/Dropout IEC 61000-4-11	>95% Dip for 0.5 Cycle 60% Dip for 5 cycles 30% Dip for 25 Cycles >95% Dip for 5 Sec	N/A. Battery Operated Device	N/A
Power Frequency 50/60 Hz Magnetic Field IEC 61000-4-8	3A/m	3A/m	Power frequency and magnetic fields should be that of a typical commercial or hospital environment.

Note: ESD compliance is dependent on a single layer of heat shrink applied over the insert earphone connector. If not installed, connector shell should not be touched during operation.

## Emissions

The Integrity System is intended for use in the electromagnetic environment specified below. The customer or user of the Integrity should ensure that it is used in such an environment.

### All equipment and systems

Emission Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1	The Integrity uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The Integrity is suitable for use in all establishments, including domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purpose.
Harmonics IEC 61000-3-2	Class N/A	
Flicker IEC 61000-3-3	N/A	

## Radio Transmissions

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment **off** and **on**, the user is encouraged to try to correct the interference by one of more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules and RSS 210 of Industry Canada. Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation of this device.

CAUTION: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

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
# MAINTENANCE INSTRUCTIONS

This product requires little maintenance. It is still recommended that a trained technician check the product every 12 months.  
For the sanitizing/disinfecting wipes, please refer to the manufacturer's cautions and instructions for use, storage, and disposal.

## AMPLITRODE® AND ELECTRODE CLIPS

The Amplitrode® and the electrode clips should be wiped with disinfecting wipes after each use.

## ABR TRANSDUCER AND ELECTRODES

### ER-3A-800

The ear tips of the ER-3A-800 Insert Earphones, and the snap electrode pads are both disposable and should be properly disposed of after each use.

### B71W Bone Conductor

The B71W and the headband should be disinfected with sanitizing wipes after each use.

## CLEANING THE OAE PROBE

For sanitary reasons and to maintain acoustic performance, the Probe must be kept clean of any visible particles, such as earwax, dead skin, and any other debris.

1. During cleaning of the Probe, hold the Probe tip down to avoid pushing the debris deeper into the sound canals of the Probe.
2. Brush off any visible debris from the Probe surface with the mini-brush.
3. Remove remaining particles from the Probe ditches with the all-purpose pick.
4. Remove particles from the sound canals and vent with the wax removal tool. The vent (a tiny hole above the cable at the back of the Probe) is to ensure the release excessive pressure from the occluded ear canal.

## DISINFECTING THE OAE PROBE

1. The Probe should be wiped with disinfecting wipes after each use. Hold the Probe tip down to avoid any ear wax or sanitizing liquid from entering the microphone sound canals.
2. Gently clean the Probe with a sanitizing wipe, and dispose of the wipe after use. Excessive force is not needed.
3. Insert the Probe tightly into the tapered hole of the Probe Holder.
4. Wipe off any excess liquid, and allow the Probe to dry for about an hour.
5. Perform a Probe Cavity Check. If a problem is detected, contact Customer Support. The problem may be due to residual ear wax on the microphone sound canals.

## ECOCHG GOLD ELECTRODE EAR TIP

Properly disposed of the gold electrode (TIProbe) ear tips after each use.

## CARING FOR THE VIVOLINK

The VivoLink may be wiped with disinfecting wipes after each use.

## CARING FOR THE COMPUTER

See the computer manufacturer's User Manual.

## DISPOSAL OF ELECTRONIC EQUIPMENT (WEEE)

The equipment shall not be disposed of as unsorted municipal waste and shall be collected separately as electrical and electronic equipment, as applicable, as specified by Directive 2012/19/EU of The European Parliament and The Council of the European Union on Waste Electrical and Electronic Equipment (WEEE).

## SYSTEM COMPONENTS

Item (Common to all Configurations)	REF/Model/Part #	Qty
System Carrying Case	VIV-BAG-1	1
Computer (with Integrity Software loaded)	M500-C	1
USB Computer Mouse	Mouse-Optical	1
Main Unit (VivoLink)	V500	1
Battery Pack	100020	2
Battery Pack Cradle	100410	1
Battery Charger	100351	1
Operator Quick Reference Guide	D-11049-2-es (este documento)	1
Super Sani Cloth wipes	H04082 (Manufacturer: PDI)	10
General Use Strap	V51-GS	1
Pediatric Shoulder Strap (2 pcs)	V51-PS	1
Item (If ABR and/or ASSR and/or 40 Hz ERP Ordered)	REF/Model/Part #	Qty
Amplitrode (in Mini Pouch)	A81 or A82	1
Insert Earphone	ER-3A-800	1
Bone Conductor	B71W	1
Ear Phone Adapter Set	ER3-06 (ABR)	1
"Snap" Skin Electrodes	Ambu <sup>®</sup> Neuroline 720 (72000-S) or Vermed <sup>®</sup> NeuroPlus (A10040 & A10041)	2 pkgs
Foam Eartips 13 mm	ER3-14A (Manufacturer: Etymotic Research)	1 pkg
Foam Eartips 10 mm	ER3-14B (Manufacturer: Etymotic Research)	1 pkg
Silicone Eartips 3.5 mm	ER3-14D2 (Manufacturer: Etymotic Research)	1 pkg
Silicone Eartips 4.0 mm	ER3-14E2 (Manufacturer: Etymotic Research)	1 pkg
Ear Tip Set (6 Each: KR006, KR036, MF-003)	ER3-V3-P (Manufacturer: Etymotic Research)	1 pkg
PDI Electrode Prep Pads	B59800 (Manufacturer: PDI)	24
Nuprep <sup>®</sup> Skin Prep Gel	10-61 (Manufacturer: Weaver and Company)	1 tube
EP Headphones (if Ordered)	H-800	1
Item (If OAE Ordered)	REF/Model/Part #	Qty
UG OAE Probe	P81-UG	1
UG OAE Probe Test Cavity	A80-14-1-0	1
UG OAE Probe Ear Tip Set	BPT3-4 (Manufacturer: Otodynamics)	1
GP OAE Probe	P81-GP	1
GP OAE Probe Test Cavity	100399	2
GP OAE Probe Ear Tip Set	TIPSET-313 (Manufacturer: Grason & Associates)	1
GP OAE Probe Cleaning Kit	HAL2573	1
Item (If ECochG Ordered)	REF/Model/Part #	Qty
ECochG Electrode Cable	ER3-60	1
Spectra 360 Electrode Gel	12-02 (Manufacturer: Parker Laboratories)	1 tube
Gold Foiled Ear Tips 13 mm	ER3-26A (Manufacturer: Etymotic Research)	1 pkg
Gold Foiled Ear Tips 10 mm	ER3-26B (Manufacturer: Etymotic Research)	1 pkg

Additional/replacement parts may be ordered in any of the following ways:



- Contact Vivosonic Customer Support at: 877.255.7685 (Canada & US) or +1-416.231.9997 (International).
- Order online at <http://www.vivosonic.com/order-accessories-howtobuy/>
- Order through their local distributor.

Please contact Vivosonic, or your local distributor, for a detailed price list.

## SPECIFICATIONS

<b>V500</b>	
Input Power	Internally powered by Vivosonic rechargeable battery pack (4 AA NiMH Batteries @ 1.2 VDC each)
Applied Part Type	BF 
Medical Device Classification	II: as per Rule 10 of Schedule 1 of SOR/98-282 - The Canadian Medical Devices Regulations IIa: as per Rule 10 of Annex IX of Council Directive 93/42/EEC Concerning Medical Devices (MDD)
Software Safety Classification	Class A – No injury or damage to health is possible (IEC 62304:2006)
Wireless Communication Type	Bluetooth: 2.402 to 2.480 MHz, Hopping Contains transmitter module FCC ID: ED9LMX9838, IC: 1520A-LMX9838 Wireless Signal Strength: 0.66 mW EIRP Emission Designation: 2M00G1D Duty Cycle: 77.02% Wireless Communication Range: 30 feet (10 meters)
Dimensions of Main Unit	L 7.1" (18 cm) x W 3.6" (9.1 cm) x H 1.2" (3.2 cm)
Weight of Main Unit	0.8 lb. (363 g) with battery pack
UMDNS Code	10-228
GMDN Category	04 Electro Mechanical Medical Devices
<b>100410 Battery Pack Cradle connected to 100351 Battery Charger</b>	
Classification	Class II
Input Power	100 - 240 VAC 50 - 60 Hz max 0.35 A

## SYMBOLS USED ON THE INSTRUMENT

Label Symbol	Description
	Read accompanying documentation
	Do not re-use. Single use only.
	Non-ionizing electromagnetic radiation (Device is capable of wireless communication)
	Electrode Contact Indicator
	Toggle ON/Standby
	Wireless communication established
+	Parking Snap for the black (+) Amplitude <sup>®</sup> clip.
1	Parking Snap for either the blue (Channel 1) Amplitude <sup>®</sup> clip or the white (-) Amplitude <sup>®</sup> clip.
2	Parking Snap for either the red (Channel 2) Amplitude <sup>®</sup> clip or the white (-) Amplitude <sup>®</sup> clip.
	Parking Snap for the green (Ground) Amplitude <sup>®</sup> clip.
<b>REF:</b>	Reference or Model Number
<b>SN:</b>	Serial Number
	Date of Manufacture
	Manufacturer
	Authorized representative in the European Community
	Type BF equipment
	Caution. Fragile. Sensitive to mechanical shock. (Usually found on packaging.)
	Keep Dry
	Input voltage
V	Voltage
	Direct Current (DC) Voltage
	Cell (battery)
AA	Cell size
NiMH	Nickel Metal Hydride battery (NiMH)
	CE Mark with notified body number
	TÜV SÜD certification
	FCC mark
	The National Institute of Metrology, Standardization and Industrial Quality (Inmetro) UL do Brasil (UL of Brazil)
	Agência Nacional de Telecomunicações Brazilian telecommunications regulator + our registration number
	The equipment shall not be disposed of as unsorted municipal waste and shall be collected separately as electrical and electronic equipment, as applicable, as specified by Directive 2012/19/EU of The European Parliament and The Council of the European Union on Waste Electrical and Electronic Equipment (WEEE).

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© 2015 Esta guía de referencia es propiedad de: